
# FORA *TN'G* BP

Blood Pressure Monitoring System

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**Owner's Manual**  
*Manual del Propietario*

# Dear FORA TN'G BP Blood Pressure Monitoring System Owner:

Thank you for purchasing the **FORA TN'G BP** Blood Pressure Monitoring System. This manual provides important information to help you operate this system correctly. Before using this product, please read the following contents thoroughly and carefully.

Managing your blood pressure is simple with this device. It's small enough to take with you anywhere, and auto-uploads all of your readings to the iFORA smartphone App. In addition, this device will help you and your healthcare professional, to keep your blood pressure under control and adjust your treatments plans as needed.

If you have other questions regarding this product, please contact the place of purchase.

# **IMPORTANT SAFETY PRECAUTIONS**

## **READ BEFORE USE**

1. Use this device **ONLY** for the intended use described in this manual.
2. Do **NOT** use the device if it is not working properly or damaged.
3. Do **NOT** use under any circumstances on newborns.
4. This device does **NOT** serve as a cure for any symptoms or diseases.  
The data measured are for reference only. Always consult your doctor to have the results interpreted.
5. Keep the equipment away from hot surfaces.
6. Do **NOT** apply the cuff to areas other than the place directed.
7. Use of this instrument in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets etc.) may cause damaging static discharges. This may cause erroneous results.
8. Do not use this instrument in close proximity to sources of strong electromagnetic radiation, as these may interfere with the accurate operation.

**KEEP THESE INSTRUCTIONS IN A SAFE PLACE**

# TABLE OF CONTENTS

<b>BEFORE YOU BEGIN</b>	<b>06</b>
Intended Use	06
Test Principle	06
Contents of System	06
Meter Overview	07
Display Screen	08
<b>HOW TO INSTALL AND UPDATE THE iFORA BP</b>	<b>09</b>
<b>BLUETOOTH PAIRING</b>	<b>10</b>
<b>IMPORTANT INFORMATION ABOUT PERFORMING A BLOOD PRESSURE TEST</b>	<b>14</b>
Before Measurement	14
Apply the Cuff Properly	14
Proper Measurement Position	16
<b>HOW TO TEST YOUR BLOOD PRESSURE</b>	<b>17</b>
Taking Measurement by the Device Itself	17
Taking Measurement with the iFORA BP App	19
<b>VIEWING PREVIOUS RESULTS</b>	<b>22</b>
Reviewing Data Record	22
Reviewing Day Average Results	24
Blood Pressure Classification	26
<b>MAINTENANCE</b>	<b>27</b>
Battery	27
Caring for Your Meter	28
<b>DETAILED INFORMATION</b>	<b>29</b>
<b>SYSTEM TROUBLESHOOTING</b>	<b>30</b>
Error Messages on Meter	30
Troubleshooting	31
<b>SPECIFICATIONS</b>	<b>32</b>
<b>FEDERAL COMMUNICATIONS COMMISSION (FCC) STATEMENT</b>	<b>33</b>

# BEFORE YOU BEGIN

## Intended Use

**FORA TN'G BP** Blood Pressure Monitoring System is a system designed to measure the systolic and diastolic blood pressure and pulse rate by using a non-invasive technique in which an inflatable cuff is wrapped on the upper arm. It is intended for individuals age 18 and up, for use in a home or clinical setting. The device is not to be used for the diagnosis of hypertension or for testing on newborns.


## Test Principle

Blood pressure is measured non-invasively at the arm based on oscillometric method.

This device is **NOT** able to take measurements in the presence of common arrhythmia, such as atrial or ventricular premature beats or atrial fibrillation. It may produce reading error.


## Contents of System

Your new **FORA TN'G BP** Blood Pressure Monitoring System includes:

- 1 Meter with Cuff
- 2 Owner's Manual
- 3 Quick Start User Guide / Warranty Card
- 4 4 x AAA Batteries
- 5 Protective Wallet

## Meter Overview

**1 BT Indicator**

Download test results with a Bluetooth connection.

**2 Cuff (immovable)**


**3 Display Screen**

**4 Battery Compartment**

**5 ON/OFF Button**

**6 BT Pairing Button**

## Display Screen

**1 Low Battery Symbol**

**5 Units of Blood Pressure**

**2 Pulse Rate Symbol**

**6 Systolic Pressure Symbol**

**3 Diastolic Pressure Value**

**7 Diastolic Pressure Symbol**

**4 Systolic Pressure Value**

**8 Pulse Rate**

# HOW TO INSTALL AND UPDATE THE iFORA BP

You can use the **FORA TN'G BP** Blood Pressure Monitoring System by itself or with the iFORA BP App.

iFORA BP is designed to assist you in easily monitoring, recording, and tracking your blood pressure. The app is simple to use and helps you to better understand your current condition as well as control your hypertension.

## System Requirement

iOS version 6 or above

Android version 4.3 or above

## Installing

iFORA BP is available on the App Store/Google Play for free download. Connect your mobile device to the internet to download the iFORA BP and then install it on your mobile device.

## Updating


When new features are added to the iFORA BP App, the App Store/Google Play will notify you automatically. Please perform the updates accordingly to provide you with the latest features that we have developed in your best interest.


# BLUETOOTH PAIRING

You can transmit your data from the meter to your devices (e.g. smart phone, tablet, PC...) via Bluetooth. Please contact your local customer service or place of purchase for assistance. Please note that you must complete the pairing between meter and Bluetooth receiver before transmitting data.

## Pairing with Your Mobile Device

1. Turn on the **Bluetooth** function on your mobile device.
2. Press the **Bluetooth button** on the meter. "PCL" will be displayed on the meter screen during the pairing process.
3. On your mobile device, launch the iFORA BP app. Select the **Setting** tab on the bottom right of the screen.


4. Tap the **Pair** button to activate the Bluetooth function.

5. Tap the **Search** button to search for available Bluetooth Smart meters. Wait for a few seconds, a pop-up window shows up to display the search results.

6. Select “TNG BP” to start pairing.

7. When pairing is complete, “TNG BP” will be displayed in the **Current Meters** field on your device.

8. Tap  to pop up the message window. Select **Yes** to save the setting.

## Bluetooth indicator on the meter

Bluetooth Indicator	Status
Flashing Blue	The Bluetooth function is enabled and waiting for connection.
Solid Blue	The Bluetooth connection is established.

### **WARNING:**

- While the meter is in transmission mode, it will be unable to perform a test.
- Make sure your device with iOS (6 or above) or Android System (4.3 or above) has the Bluetooth function turned on and the meter is within the receiving range before transmitting the data.
- A compatibility issue between your mobile device and the meter might occur due to the way Bluetooth functionality is implemented by the various mobile device manufacturers.


# IMPORTANT INFORMATION ABOUT PERFORMING A BLOOD PRESSURE TEST

## Before Measurement

- Avoid caffeine, tea, alcohol and tobacco for at least 30 minutes before measurement.
- Wait for 30 minutes after exercising or bathing before measurement.
- Sit or lie down for at least 10 minutes before measurement.
- Do not measure when feeling anxious or tense.
- Take a 5-10 minute break between measurements. This break can be longer if necessary, depending on your physical condition
- Keep the records for your doctor as reference.
- Blood pressure naturally varies between each arm. Always measure your blood pressure on the same arm.

## Apply the Cuff Properly

- 1.** Stretch your left (right) arm in front of you with your palm facing up. Slide and place the cuff onto your arm, the arrow should point toward the upper arm.
- 2.** Wrap and tighten the cuff above your elbow. The red line on the edge of the cuff should be approximately **0.8" to 1.2" (2 cm to 3 cm)** above your elbow.

**3.** Leave a little free space between the arm and the cuff; you should be able to fit two fingers between them. Clothing must not restrict the arm. Remove all clothing covering or constricting the measurement arm.


**4.** Press the hook material firmly against the pile material. The top and bottom edges of the cuff should be tightened evenly around your upper arm.

**NOTE:**

The cuff is fixed on the blood pressure monitor. Do **NOT** remove the cuff or change the size of cuff.

## Proper Measurement Position

1. Sit down for at least 10 minutes before measuring.
2. Place your elbow on a flat surface. Relax your hand with the palm facing up.
3. Make sure the cuff is about the same height as the location of your heart.
4. Remain still and do not talk or move during the measurement.

### **WARNING:**

If the cuff is relatively lower (higher) than the heart, the obtained blood pressure value could be higher (lower) than the actual value. A 15 cm difference in height may result in an error around 10 mmHg.

# HOW TO TEST YOUR BLOOD PRESSURE

## Taking Measurement by the Device Itself

**Always apply the pressure cuff before turning on the meter.**

**1. Tap ON/OFF Button.** All the LCD symbols will appear with a long “beeping” sound. Then the cuff will begin to inflate automatically.

**2. Pump starts to inflate.** The pressure in cuff increases to working pressure.

**3. The heart symbol “♥” will show** silently when pulse is detected during the measurement. When signal of pulse becomes strong, “♥” icon starts to flash.

### NOTE:

When the device detects irregular heart beat, the IRB symbol is shown as “”.

**This is for your reference only; if the problem still persists or if you feel uncomfortable, please consult your doctor.**

4. After the measurement, the meter displays the systolic pressure, diastolic pressure and pulse rate.

5. Tap **ON/OFF Button** to switch off. The meter will switch off automatically after being idle for 3 minutes.


#### NOTE:

- If you tap **ON/OFF Button** during measurement, the meter will be turned off.
- When the “Syn” icon appears on the screen, it indicates you have not paired the meter to your mobile device. You need to synchronize with your mobile device in order to adjust the date and time settings. In addition, this measurement will not be stored in meter memory.


## Taking Measurement with the iFORA BP App

All the screenshots displayed below come from the display screen of the Bluetooth receiver installed with iFORA BP App, not from the display screen of the meter.

### Single Measurement

1. Turn on the Bluetooth switch on your mobile device and launch iFORA BP App.
2. Apply the pressure cuff, and then tap **ON/OFF Button** on the iFORA BP screen. All the LCD symbols will appear with a long “beeping” sound on the device monitor. The cuff will begin to inflate automatically until it reaches the working pressure.
3. The pulse rate is detected during the measurement.
4. After the measurement, the meter displays testing results of systolic pressure, diastolic pressure and pulse rate.
5. Tap Pulse rate, a predictive analysis page is shown on the screen.
6. The meter will switch off automatically after being idle for 3 minutes.


## Average Measurement

This function automatically performs three (3) consecutive blood pressure measurements and displays the average result in the end. This function is only available when your meter is connected to your mobile device.

**1.** Press the **Bluetooth button**. The “PCL” icon will appear on the meter display, indicating that the device is ready to transmit data.

**2.** Launch iFORA BP App.

**3.** Enable **Average Mode** and tap **ON/OFF Button** on the iFORA BP screen. Then the cuff will begin to inflate automatically.


**4.** After the first measurement is complete, the device and iFORA BP app will start counting down before the second measurement begins. The device will take three (3) measurements consecutively with an interval of 20 seconds.

### NOTE:

When the difference between the first and second systolic pressure is more than 15 mmHg, the time interval for third measurement will be 40 seconds.

5. After the measurement, the average result is displayed on the screen with an “**AVG**” tag.

6. Tap Pulse rate, a predictive analysis page is shown on the screen.

7. The meter will switch off automatically after being idle for 3 minutes.

**NOTE:**

- If you tap **ON/OFF Button** on the iFORA BP screen during measurement, the meter will be turned off.
- Always apply the pressure cuff before turning on the meter.


# VIEWING PREVIOUS RESULTS

## Reviewing Data Record

Your meter stores the 200 most recent blood pressure test results along with respective dates and times in the memory. To recall the memory, you have to connect your **FORA TN'G BP** Blood Pressure Monitoring System to your mobile device. Once connected, the measured data is automatically imported to iFORA BP App.

Tap “**Data**” on the main menu bar to review the records. The first reading you see is the latest blood pressure result.

### Data Record Screen

#### NOTE:

- Swipe horizontally to switch between Morning & All measurements.
- Tap the “Blood pressure classification” for more details (see Blood pressure classification).

## Data Record Screen (rotate by 90°)

Rotate by 90° to see your results sorted by Date / Time.

Measurement mode

Irregular heart beat

Date/Time	Mode	Blood Pressure	IRB	Pulse	Note
Apr 24, 2015 2:13 PM		104 / 62 mmHg		73 bpm	
Apr 20, 2015 6:43 PM		119 / 81 mmHg		69 bpm	
Apr 16, 2015 11:21 AM		89 / 54 mmHg		85 bpm	
Apr 16, 2015 11:19 AM	 AVG	91 / 54 mmHg		76 bpm	
Apr 16, 2015 11:14 AM		89 / 53 mmHg		77 bpm	
Apr 13, 2015 6:10 PM		111 / 88 mmHg		73 bpm	

Systolic / Diastolic blood pressure

Sunny symbol indicates Morning measurement

### NOTE:


Scroll vertically to view history records.

## Reviewing Day Average Results

iFORA BP offers you 7-, 14-, and 30- day average results to assist you to monitor your blood pressure.

Tap “**Analysis**” on the main menu bar to enter Day Average Results.

### Day Average Screen


#### NOTE:

- Swipe horizontally to switch between Morning & All measurements.
- Tap the “Blood pressure classification” for more details (see Blood pressure classification).

## Day Average Screen (rotate by 90°)

Rotate by 90° to see your results in trend graph.

### NOTE:

- Scroll horizontally to view older results.
- The function of viewing earlier results is featured in the Bluetooth receiver installed with the iFORA BP.
- The meter itself only displays the latest result after each measurement.

## Blood Pressure Classification

Red

Orange

Green



Light Green

Yellow

Under Data Record screen or Day Average screen, you can tap on the Blood pressure classification to enter the single or total result classification.

Here, your blood pressure results are classified by colors based on your Systolic and Diastolic readings, with Red for Stage 2 Hypertension, Orange for Stage 1 Hypertension, Green for Prehypertension, Light Green for Normal and Yellow for too low. The blood pressure classification is based on JNC VII Express (2003).

### NOTE:

- The information here is only displayed on the screen of Bluetooth receiver installed with the iFORA BP App.
- The classification is for reference only. Always consult your doctor to have the results interpreted.

# MAINTENANCE

## Battery

Your meter comes with four (4) 1.5 V AAA size alkaline batteries.

### Low Battery Signal

The meter will display either of the two messages below to alert you when the meter power is getting low.

**1. The  symbol appears** along with display messages: The meter is functional and the result remains accurate, but it is time to change the batteries.

**2. The  symbol appears with E-b:** The power is not enough to do a test. You must change the batteries immediately.

### Replacing the Battery

**To replace the batteries, make sure the meter is turned off.**

- 1.** Press the edge of the battery cover and lift it up to remove.
- 2.** Remove the used batteries.
- 3.** Replace with four 1.5 V AAA size alkaline batteries in accordance with the mark +/- indication in the battery compartment.
- 4.** Close the battery cover.

**NOTE:**

- As with all small batteries, these batteries should be kept away from small children. If swallowed, promptly seek medical assistance.
- Batteries might leak chemicals if unused for a long time. Remove the batteries if you are not going to use the device for an extended period (i.e., 3 months or more).
- Properly dispose the batteries according to your local environmental regulations.
- Replacing the batteries does not affect the test results stored in memory. However, if you did not synchronize with your mobile device after the battery replacement, the measured data will not be stored in the memory.

## Caring for Your Meter

To avoid the meter from attracting dirt, dust or other contaminants, wash and dry your hands thoroughly before use.

### Cleaning

- To clean the meter exterior, wipe it with a cloth moistened with tap water or a mild cleaning agent, then dry the device with a soft dry cloth. Do NOT flush with water.
- Do NOT use organic solvents to clean the meter.
- Do NOT wash or iron the pressure cuff.

### Meter Storage

- Storage condition:  $-25^{\circ}\text{C}$  to  $70^{\circ}\text{C}$  ( $-13^{\circ}\text{F}$  to  $158^{\circ}\text{F}$ ), 10% to 95% relative humidity.
- Always store or transport the meter in its original storage case.
- Avoid dropping or heavy impact.
- Avoid direct sunlight and high humidity.

# DETAILED INFORMATION

## Reference Values

Clinical studies show that adult diabetes is often accompanied by elevated blood pressure. People with diabetes can reduce their heart risk by managing their blood pressure along with diabetes treatment\*1.

Monitoring your routine blood pressure trend helps you to know your body condition. Human blood pressure naturally increases after reaching middle age. This symptom is a result of continuous ageing of the blood vessels. Further causes include obesity, lack of exercise and cholesterol (LDL) adhering to the blood vessels. Rising blood pressure accelerates hardening of the arteries, and the body becomes more susceptible to apoplexy and coronary infarction. The recommended blood pressure range is as below:

Classification	Systolic Pressure (mmHg)	Diastolic Pressure (mmHg)
Hypotension *2	≤ 90	≤ 60
Normal *3	< 120	< 80
Pre-hypertension *3	120–139	80–89
Stage 1 Hypertension *3	140–159	90–99
Stage 2 Hypertension *3	≥ 160	≥ 100

\*1: American Diabetes Association: The Diabetes-Heart Disease Link Surveying Attitudes, Knowledge and Risk (2002).

\*2: National Heart, Lung, and Blood Institute, Diseases and Conditions.

\*3: The Seventh Report of the Joint National Committee on Prevention, Detection, Evaluation, and Treatment of High Blood Pressure. NIH Publication. 2003. No. 03-5233.

# SYSTEM TROUBLESHOOTING

If you follow the recommended action but the problem persists, or error messages other than the ones below appear, please call your local customer service. **Do not attempt to repair by yourself and never try to disassemble the meter under any circumstances.**

## Error Messages on Meter

MESSAGE	CAUSE	WHAT TO DO
	Inflation or pressure error.	Please contact local customer service for help.
	Blood pressure measurement error.	Refit cuff tightly and correctly. Relax and repeat the measurement. If error still remains, contact local customer service for help.
 	Problems with the meter.	Review the instructions and repeat the test. If the meter still does not work, please contact the local customer service for help.
	Battery is too low.	Repeat with new batteries.

### NOTE:

The information here is only displayed on the meter screen.

## Troubleshooting

### 1. If no display appears after tapping **ON/OFF Button**:

POSSIBLE CAUSE	WHAT TO DO
Batteries exhausted.	Replace the batteries.
Batteries incorrectly installed or absent.	Check that the batteries are correctly installed.

### 2. If the heart rate is higher/lower than user's average:

POSSIBLE CAUSE	WHAT TO DO
Movement during measurement.	Repeat measurement.
Measurement taken just after exercise.	Rest for at least 30 minutes before repeating measurement.

### 3. If the result is higher/lower than user's average measurement:

POSSIBLE CAUSE	WHAT TO DO
May not be in correct position while measuring.	Blood pressure naturally varies from time to time.
Blood pressure naturally varies from time to time.	Keep in mind for next measurement.

### 4. If the cuff inflates again while measuring:

POSSIBLE CAUSE	WHAT TO DO
Cuff is not fastened.	Cuff is not fastened.
If user's blood pressure is higher than the pressure the device has inflated, the device will automatically increase the pressure and start to inflate again. Stay relaxed and wait for the measurement.	

# SPECIFICATIONS

## ● System Performance

**Power Source:** Four 1.5 V AAA alkaline batteries

**Size of Meter w/o Cuff:**

145 (L) x 64.5 (W) x 29.6 (H) mm, 240 g (without batteries)

**Power Saving:** Automatic power off if system idle for 3 minutes

**Waterproof Grade:** IP22

**The Expected Service Life:** 2 years

**System Operating Condition:**

41 °F to 104 °F (5 °C to 40 °C), 15% to 93% RH, 700 hPa to 1060 hPa

**Meter Storage Condition:**

-13 °F to 158 °F (-25 °C to 70 °C), 10% to 95% RH

## ● Blood Pressure Measurement Performance

**Pressure Range:** 0 - 300 mmHg

**Systolic Measurement Range:** 50 mmHg -250 mmHg

**Diastolic Measurement Range:** 30 mmHg -180 mmHg

**Pulse Rate Measurement Range:** 40 -199 beats / minute

**Maximum Inflation Pressure:** 280 mmHg

**Accuracy of Pressure:**  $\pm 3$  mmHg or  $\pm 2\%$  of reading

**Accuracy of Pulse Rate:**  $\pm 4\%$  of reading

**Measurement Unit:** mmHg

This device has been tested to meet the electrical and safety requirements of: IEC 60601-1, and IEC 60601-1-2

**Reference to Standards:**

- IEC60601-1 General requirement for safety
- IEC60601-1-2 Requirements for EMC

**IP22:**

- Object size protected against:  $>12.5$  mm (Fingers or similar objects)
- Protected against: Dripping water when tilted up to  $15^\circ$ .

# FEDERAL COMMUNICATIONS COMMISSION (FCC) STATEMENT

15.21

**You are cautioned that changes or modifications not expressly approved by the part responsible for compliance could void the user's authority to operate the equipment.**

15.105(b)

## **Federal Communications Commission (FCC) Statement**

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

**This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:**

- 1) This device may not cause harmful interference and
- 2) This device must accept any interference received, including interference that may cause undesired operation of the device.

## FCC RF Radiation Exposure Statement:

- 1) This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.
- 2) This equipment complies with FCC RF radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with a minimum distance of 20 centimeters between the radiator radiation source and your body.

<b>Guidance and manufacturer's declaration-electromagnetic emissions</b>		
The <b>FORA TN'G BP</b> Blood Pressure Monitoring System is intended for use in the electromagnetic environment specified below. The customer or the user of the <b>FORA TN'G BP</b> Blood Pressure Monitoring System should assure that it is used in such an environment.		
<b>Emission test</b>	<b>Compliance</b>	<b>Electromagnetic environment-guidance</b>
RF emissions CISPR 11	Group 1	The <b>FORA TN'G BP</b> Blood Pressure Monitoring System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The <b>FORA TN'G BP</b> Blood Pressure Monitoring System is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Not applicable	

### Guidance and manufacturer's declaration-electromagnetic immunity

The **FORA TN'G BP** Blood Pressure Monitoring System is intended for use in the electromagnetic environment specified below.

The customer or the user of the **FORA TN'G BP** Blood Pressure Monitoring System should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+ 6 kV contact + 8 kV air	+ 6 kV contact + 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC 61000-4-4	+ 2kV for power supply lines + 1kV for input/ output Lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+ 1kV line(s) to line(s) + 2kV line(s) to earth		
Voltage Dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0,5 cycle 40% UT(60% dip in UT) for 5 cycles 70% UT(30% dip in UT) for 25 cycles <5% UT(>95% dip in UT) for 5 s	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the <b>FORA TN'G BP</b> Blood Pressure Monitoring System requires continued operation during power mains interruptions, it is recommended that the <b>FORA TN'G BP</b> Blood Pressure Monitoring System be powered from an uninterruptible power supply or a battery.

Power frequency (50/60 HZ) magnetic field IEC 61000-4-8	3 A/m	3 A/m	The <b>FORA TN'G BP</b> Blood Pressure Monitoring System power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: UT is the a.c. mains voltage prior to application of the test level.			

<b>Guidance and manufacturer's declaration-electromagnetic immunity</b>			
<p>The <b>FORA TN'G BP</b> Blood Pressure Monitoring System is intended for use in the electromagnetic environment specified below.</p> <p>The customer or the user of the <b>FORA TN'G BP</b> Blood Pressure Monitoring System should assure that is used in such and environment.</p>			
<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment-guidance</b>
Conducted RF IEC 61000-4-6	3 Vrms 150 KHz to 80 MHz	Not applicable	<b>Portable and mobile RF communications equipment should be used no closer to any part of the FORA TN'G BP</b> Blood Pressure Monitoring System including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	Recommended separation distance: $d = 1,2 \sqrt{P}$ $d = 1,2 \sqrt{P}$ 80MHz to 800 MHz $d = 2,3 \sqrt{P}$ 800MHz to 2,5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).

			Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup> , should be less than the compliance level in each frequency range <sup>b</sup> . Interference may occur in the vicinity of equipment marked with the following symbol: 
--	--	--	---

NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 1: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the **FORA TN'G BP** Blood Pressure Monitoring System is used exceeds the applicable RF compliance level above, the **FORA TN'G BP** Blood Pressure Monitoring System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the **FORA TN'G BP** Blood Pressure Monitoring System.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**Recommended separation distance between portable and mobile RF communications equipment and the FORA TN'G BP Blood Pressure Monitoring System**

The **FORA TN'G BP** Blood Pressure Monitoring System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the **FORA TN'G BP** Blood Pressure Monitoring System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the **FORA TN'G BP** Blood Pressure Monitoring System as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz $d = 1,2\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3\sqrt{P}$
0,01	N/A	0,12	0,23
0,1	N/A	0,38	0,73
1	N/A	1,2	2,3
10	N/A	3,8	7,3
100	N/A	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $p$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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**Foracare, Inc.**