

尺寸: 440*290mm 材质: 80g书纸, 单黑折页 (风琴折再对折)


——— 切线
- - - - - 折痕线

Andesfit health

Instruction manual

Infrared Forehead Thermometer

Model:UFR103



Introduction

▲ Your new Infrared forehead thermometer uses advanced infrared (IR) technology to measure temperature instantly and accurately on the forehead or object.

▲ Easy to use and less measurement time
This thermometer does not need to contact body or object to ensure the safety and hygiene. Its ergonomic design makes this thermometer be simple and very easy to check the temperature. It only takes 1 second to take measurement and reading.

▲ Body and object modes

This thermometer supports to measure temperature of body and object. Measuring range of object mode is 0.0 to 100.0°C (32.0°F-199.9°F). That means except to take body temperature, it also supports to take the surface temperature of object as below:

Surface temperature of milk in the bottle;

Surface temperature of baby bath;

Environment temperature.

▲ Alarm sounds

When body temperature is over 37.5°C, an alarm sounds for 10 times.

▲ Memory function

It stores up to 20 sets recent measurement data.

▲ Auto power off

Automatically power off if left idle for Approx 60±5s.

▲ Measurement

Measurement time interval in 1s and the measurement distance is within 1-3cm. Please read the manual carefully before you use the unit, and keep for future reference.

▲ Intended use:

The Infrared Thermometer is intended for the intermittent measurement and monitoring of human body temperature from forehead or object. The device is indicated for use by people of all ages at homecare and in hospital.

Safety Information

To assure the correct use of the product, basic safety measures should always be followed including the warning and the caution listed in the instruction manual:

Symbol descriptions

The following symbols may appear in this manual, on the label, on the device, or on its accessories. Some of the symbols represent standards and compliances associated with the device and its use.

	WARNING: This alert identifies hazards that may cause serious personal injury or death.
	CAUTION: This alert identifies hazards that may cause minor personal injury, product damage, or property damage.
	Type BF applied part
	Manufacturer
	Specifies serial number
	Authorized Representative in the European Community
	CE Mark: conforms to essential requirements of the Medical Device Directive 93/42/EEC.
	DISPOSAL: Do not dispose this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary.
	Direct current
	Follow instructions for use
	CAUTION: Consult accompanying documents

▲ This thermometer is not intended to substitute for a consultation with your physician. The Forehead scan temperature serves as a reference only. It cannot be a judgment on fever.

▲ Basic safety precautions should always be observed, especially when the thermometer is used on or near children and disabled persons.

▲ Please place the device unreachable by young.

▲ Avoid direct sunlight.

▲ Do not touch the lens.

▲ No modification of this device is allowed.

▲ The swallowing of small parts like packing bag, battery, battery cover and so on may cause the suffocation.

▲ Please do not use a dilution agent, alcohol or petrol to clean the unit. Please treat it gently and prevent the falling from a high place.

▲ Please do not immerse it in liquid.

▲ Never leave battery in the battery compartment for a long time without use, as they may leak and cause damage to the unit.

▲ Please take off the battery if you do not intend to use within 3 months. Replace with new batteries if the unit display a low battery symbol.

▲ Do not mix the old and new batteries together.

▲ Do not use during the transportation.

▲ WARNING:

Do not dispose of electrical appliances as unsorted municipal waste. Use separate collection facilities. Contact your local government for information regarding the collection systems available. If electrical appliances are disposed of in landfills or "burning," hazardous substances can leak into the groundwater and get into the food chain, damaging your health and well-being.

Classification


- Internally powered equipment;
- Type BF applied part;
- Protection against ingress of water or Particulate matter: IP21;
- Not category AP (APG) equipment;
- Mode of operation: Continuous operation;

▲ The user must check that the equipment functions safely and see that it is in proper working condition before being used.


Product Structure

-Display

- Object mode
- Forehead temperature mode
- Ambient temperature mode
- Temperature unit: °C/°F
- Temperature reading
- Memory symbol
- Battery symbol
- Horn symbol
- Number of memory groups
- Bluetooth symbol



-Body



Battery Installation

Remove the battery cover from the battery compartment, insert the battery.

- Remove the battery cover from the battery compartment as the arrow direction accordingly.



- Insert 2 AAA powerful batteries into the compartment and ensure each battery is in the proper direction. Positive (+) and Negative (-) are displayed on the back of battery cover.

- Close the battery cover.

Low battery and replacement

When power on, the low battery symbol will display once the unit start to work, and you must replace with new batteries, otherwise the unit can't work.

Battery type and replacement

Please use 2pcs AAA identical 1.5V alkaline batteries.

Do not use the batteries beyond their expiry date.

Please remove the batteries if you do not need to use for long time.

WARNING


Dispose of the battery in accordance with all federal, state and local laws. To avoid fire and explosion hazard, do not burn or incinerate the battery.

Setting Mode

How To Set

- Mode setting:

The display all shows for one second after turning on the On-off key, the default setting is "Body mode" after turning on. Short press the setting button once to switch to "Object mode", and then press the setting button again to switch to "Ambient temperature mode". If you press the setting button again, it will return to the "Body mode".



Note: The Body Mode is used to measure the forehead temperature, and the Object Mode is used to measure the object temperature, then the Ambient mode is used to measure the Ambient temperature.

- Unit setting and sound switch settings:

Long press the setting button for 3 seconds to enter the °C/°F and sound switch settings.

Step 1: Enter °C / °F switch: the interface displays °C or °F, press the memory key to switch the temperature unit, press the set key to confirm the selection, at the same time enter the sound setting selection mode. (Restart, the unit will default to °C).



Step 2: Enter the sound switch setting: the interface displays On or OFF, press the memory key to switch. If you select on, the horn symbol is displayed on the screen. If you select OFF, the horn symbol is not displayed on the side. Press the set button to confirm the selection. The setting is completed automatically. After setting, it will automatically return to the measure mode and wait for measurement. (Restart, the sound will default to on)



Proper Use Of The Unit

Pre-measurement

About Normal Body Temperature&Fever

The temperature in the forehead and temple area differs from the internal temperature, which is taken orally or rectally. Vasoconstriction, an effect which constricts the blood vessels and cools the skin, can occur during the early stages of a fever. In this case, the temperature measured by the Infrared thermometer may be unusually low. If the measurement therefore does not match the patient's own perception or is unusually low, re-peat the measurement every 15 minutes. As a reference, you can also measure the internal body temperature using a conventional oral or rectal thermometer.

Body temperature can vary from one individual/person to next. It also varies by location on the body and time of day. Below shows the statistical normal ranges from different sites.

Please keep in mind that temperatures measured from different sites, even at the same time, should not be directly compared. Fever indicates that the body temperature is higher than normal. This symptom may be caused by infection, over-dressing or immunization. Some people may not experience fever even when they are ill.

These include, but are not limited to, infants younger than 3 months old, persons with compromised immune systems, persons taking antibiotics, steroids or antipyretics (aspirin, ibuprofen, acetaminophen), or persons with certain chronic illnesses. Please consult your physician when you feel ill even if you do not have fever.

Table*1 Body Site Normal Temperature Range

Body Site	Normal Temperature Range
Oral	0.6°C(1°F) or more above or below 37°C(98.6°F)
Rectal/ear	0.3°C to 0.6°C (0.5°F to 1°F) higher than oral temperature
Axillary (armpit)	0.3°C to 0.6°C (0.5°F to 1°F) lower than oral temperature

Note: Body Temperature at WebMD, website: <http://firstaid.webmd.com/body-temperature>, retrieved at 2010 Jan 7.

As a Body Thermometer

Press the On-off button, after the display is fully displayed for 1 second, the measurement preparation screen is displayed. The C icon lights up. Align the unit with the center of the forehead and let the infrared thermometer probe be 1cm away from the forehead, measurement completion prompts measurement results. (Buzzer: One long beep Voice)



Note:

- If the reading is $\geq 37.5^{\circ}\text{C}$ (99.5°F) and $< 43^{\circ}\text{C}$ (109.4°F), the buzzer ten short beeps.
- After the measurement is completed, the measurement can be continued after 1 second.
- As the forehead measurement temperature is likely to be affected by sweat, oil and the surroundings, the reading shall be taken as a reference only.
- If the probe is placed at an angle close to the forehead measurement, the reading will be affected by surrounding temperature.
- Babies' skin reacts very quickly in the ambient temperature. Therefore, do not take their temperature with the non-contact thermometer during/after breastfeeding, because the skin temperature may be lower than the internal body temperature.
- The user and the infrared thermometer should be quiescent for at least 30 minutes under similar indoor conditions.

As an Object Thermometer

Press the On-off button, after the display is fully displayed for 1 second, switch to the object mode and display the measurement preparation screen. The °C icon will light up. Point the infrared thermometer at the object to be measured, and press On-off button to measure, measurement completion prompts measurement results.(Buzzer: One long beep Voice)



Note:

- As the object scan temperature is likely to be affected by the surroundings, the reading shall be taken as a reference only.
- If the probe is placed at an angle close the object scan, the reading will be affected by surrounding temperature.

As an Ambient Thermometer

Press the On-off button, after the display is fully displayed for 1 second, switch to the ambient temperature mode, display the measurement preparation screen, the °C icon lights up, press the On-off key measurement, measurement completion prompts measurement results.(Buzzer: One long beep Voice)



Memory-recall of measurements

When the device is off, press the "memory button" to enter the read memory mode, the memory icon and the first group of memory values appear. Each time you press the key in sequence, 20 groups of memories are read in sequence, and the first group of memories The value is the latest measurement reading, and so on. When reading the 20th group of memory, if you press the memory key again, the first group of readings will be re-read.



Memory-clear Of Measurements

In the shutdown state, long press the memory button for about 8 seconds, the LCD screen displays CL, press the On-off key to confirm, CL flashes 3 times to clear all memories, and the buzzer beeps 3 times continuously



Exceptional Situation

Symbol	Correction
	In Body Mode, Measured temperature is above Measuring range of 43°C / 109.4°F it is shown on the LCD.

	In Body Mode, Measured temperature is below the Measuring range of 32°C / 89.6°F , it is shown on the LCD.
	In Object Mode, Measured temperature is above Measuring range 100.0°C / 198.9°F or environmental temperature is above system Operating range 40°C / 104.0°F , it is shown on the LCD.
	In Object Mode, Measured temperature is below Measuring range 0.0°C / 32.0°F or environmental temperature is below system operating range 5°C / 41.0°F , it is shown on the LCD.
	In Ambient Mode, Measured temperature is above Measuring range 40.0°C / 104.0°F , it is shown on the LCD.
	In Ambient Mode, Measured temperature is below Measuring range 0.0°C / 32.0°F , it is shown on the LCD.
	Low battery, replace all the worn batteries with new ones.
	Thermometer system fails or affected by electric magnetic field.

▲ Please contact the distributor if you can not solve the problem, do not disassemble the unit by yourself!

Care And Maintenance

Care for the main unit

- Keep the unit in the storage case when not in use.
- Clean the unit with a soft dry cloth. Do not use any abrasive or volatile cleaners.
- Never immerse the unit or any of its' component in water.

Maintenance

Do not clean the unit with naptha, thinner or gasoline etc. Store the unit in a clean and dry location, not subject the unit to extreme hot or cold temperature, humidity and direct sunlight.

	Remove the batteries if the unit will not be used in 3 months or longer.		Do not use the unit under the influence of electromagnetic interference (nearby cellphones, microwave etc.)
--	--	--	---

Note: We will not be responsible for any quality problem if you do not care and maintain the product as instructed.

Specification

Description	Non-contact Infrared Forehead Thermometer
Display	LCD digital display
Measuring location	Forehead and object surface
Measurement range	Body mode : 32.0°C - 43.0°C (89.6°F - 109.4°F) Object mode : 0.0°C - 100.0°C (32.0°F - 199.9°F) Ambient temperature : 0.0°C - 40.0°C (32.0°F - 104.0°F)
Temperature unit	°C/°F
Display resolution	$0.1^{\circ}\text{C}/0.1^{\circ}\text{F}$
Accuracy	$\pm 0.2^{\circ}\text{C}/\pm 0.4^{\circ}\text{F}$ (within 35.0°C - 42.0°C / 95.0°F - 107.6°F)
Memory function	20 sets memory of measurement values
Buzzer function (or voice reading)	(1) Turn on the device : 1short beep (2) Measurement completed: 1long beep (3) Fever: 37.5°C or 99.5°F : 10 short beeps
Backlight	Green
Power source	2*1.5V AAA (2 pcs AAA alkaline battery)
Automatic power-off	In 60s.5s
Main unit weight	Approx.98g (batteries not included)
Main unit size	115*22mm*W30mm*H38mm
Battery life	Could be used for 300 times for normal condition

Accessories	Instruction manual
Operating environment	Body mode: 10°C - 40°C (50°F to 104°F) Object mode: 5°C - 40°C (41°F to 104°F) Relative humidity range : $\leq 85\%$ RH Atmospheric pressure range : 750hPa - 106hPa .
Storage environment	Ambient temperature range : -20°C - 50°C Relative humidity range : 15% - 95% RH Atmospheric pressure range : 750hPa - 106hPa .

Clinical Measurement Accuracy And Safety Verification:

The product has passed clinical trials, The measured results of the infrared forehead thermometer was compared with the measured results of mercury thermometers, the deviation average $\Delta t = 0.011^{\circ}\text{C}$ not exceeding 0.3°C ; the clinical repeatability of the infrared forehead thermometer $\text{SR} = 0.100^{\circ}\text{C}$, not exceeding 0.3°C . The measured results up to the laboratory standard and the clinical standard, Therefore, the deviation average and the clinical repeatability of the infrared forehead thermometer are complied with the regulatory requirement ISO 80601-2-56. The conclusions are drawn from the clinical trials, the accuracy and safety are compliedwith the regulatory requirement.

Warranty Information

▲ The unit is guaranteed to be free of defects in workmanship and materials under normal use for a period of 1 Year from the date of purchase.

▲ For repair under this warranty, Our authorized service agent must be advised of the fault within the period of the warranty. This warranty only covers parts and labor service under normal operations. Any defect resulting

from natural causes, eg, flood, hurricane etc, is not covered in this guarantee. This guarantee also does not cover damage incurred by use of the unit not in accordance with the instructions, accidental damage, or being tampered or serviced by unauthorized service agents.

▲ The flowing will be excluded from this warranty-if the thermometer has been misused, abused, or neglect in following the manual's instructions on purpose and unauthorized repair or modifications.

▲ The device requires no calibration.

▲ The device is not repairable and contains no user serviceable parts.

EMC

IEC 60611-2, 2014 ME EQUIPMENT and ME SYSTEMS identification, marking and documents for Class B product

Instructions for use

The ME EQUIPMENT or ME SYSTEM is suitable for home healthcare environments and so on.

Warning: Don't near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and affect its proper operation.

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Infrared Ear Thermometer (TE-66, TE-68, TE-82), including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

If any: a list of all cables and maximum lengths of cables (if applicable), transducers and other ACCESSORIES that are replaceable by the RESPONSIBLE ORGANIZATION and that are likely to affect compliance of the ME EQUIPMENT or ME SYSTEM with the requirements of Clause 7 (EMISIONS) and Clause 8 (IMMUNITY), ACCESSORIES may be specified either generically (eg, shielded cable, lead impedance) or specifically (e.g., by MANUFACTURER and EQUIPMENT OR TYPE REFERENCE).

If any: the performance of the ME EQUIPMENT or ME SYSTEM that was determined to be ESSENTIAL PERFORMANCE and a description of what the OPERATOR can expect if the ESSENTIAL PERFORMANCE is lost or degraded due to EM DISTURBANCES (the defined term "ESSENTIAL PERFORMANCE" need not be used).

Technical description

1. All necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the expected service life.

2. Guidance and manufacturer's declaration -electromagnetic emissions and Immunity

Guidance and manufacturer's declaration - electromagnetic emissions	
Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	N/A
Voltage fluctuations / flicker emissions IEC 61000-3-3	N/A

Guidance and manufacturer's declaration - electromagnetic immunity		
Immunity Test	IEC 60601-1-2 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Electrical fast transient/burst IEC 61000-4-4	Power supply lines: ± 2 kV Signal/buss lines: ± 1 kV	N/A
Surge IEC 61000-4-5	Line(s) to line(s): ± 1 kV Line(s) to earth: ± 2 kV, 100 kHz repetition frequency	N/A
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% 0.5 cycle AC: 25° - 30°C , 195° - 190°C DC: 225° - 270° and 315° 0% 1 cycle And 70% 25/30 cycles Single phase: at 0 0% 300 cycle	N/A
Power frequency magnetic field IEC 61000-4-8	30 A/m 50/60Hz	30 A/m 50/60Hz
Conducted RF IEC 61000-4-6	150kHz to 80MHz: 30V/m 0.1V/m in ISM and amateur radio bands 30% Am at 1kHz	N/A
Radiated RF IEC 61000-4-3	10 V/m 80 MHz - 2.7 GHz: 80 $\mu\text{V}/\text{m}$ at 1 kHz	10 V/m 80 MHz - 2.7 GHz: 80 $\mu\text{V}/\text{m}$ at 1 kHz

NOTE: UT is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - electromagnetic immunity							
Test Frequency (MHz)	Band (MHz)	Service	Modulation (V)	Modulation (V)	Distance (m)	IMMUNITY TEST LEVEL (V/m)	
Radiated RF IEC 61000-4-3 (Test specification for SURE PORT RF wireless communication equipment)	385 - 390 450 - 300 745 - 767 750 - 810 810 - 860 870 - 930 1720 - 1900 1645 - 1970 2450	TETRA 400 GSM-R 460, FRS 460 GSM 800/900, LTE Band 3, 4, 25, UMTS GSM 1900, CDMA 1900, GSM 1900, DECT, LTE Band 1, 3, 4, 25, UMTS Bluetooth, WLAN 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 18 Hz FM ± 5 kHz deviation 1 kHz sine Pulse modulation 217 Hz Pulse modulation 18 Hz Pulse modulation 217 Hz	1.8	0.3	27	
	5240 - 5100 5000 - 5785	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9	

Andesfit Health
Shenzhen Andesfit Technology Co., Ltd.
Floor 404, Building 1, 19th Street, Technology Industrial Zone, No.1, Chuanxin Road, Heshui Community, Matian Street, Guangzhou, New District, 510100, Guangdong, P.R. CHINA
Tel: 0086-20-22212108 E-mail: andesfit@andefit.com
MADE IN CHINA
Shenzhen Andesfit Technology Co., Ltd. (Shenzhen Office)
Tel: 0086-755-23111515 E-mail: andesfit@andefit.com

Andesfit Ltd
Unit 513 Lakeside 1, No.3 Science Park West Ave.,
Hong Kong Science Park, Hong Kong
Customer Service: 0012 3591380
www.andesfit.com

